

[Product Range](#) | [Polymers](#) | [Ordering + Selection](#) | [Fitting](#) | [Maintenance + Storage](#) | [Location](#)

[Sieves](#) | [Filters](#)

[Home](#) > [Product Range](#) > [Sieves](#) > [MiniSIV](#)

 [Printer Friendly Version](#)

MiniSIV

Normesh MiniSIV product features:

Tel: 0161-628 9849
Fax: 0161-627 5732

E-mail: [SALES](#)
[PURCHASING](#)
[ACCOUNTS](#)

Normesh Limited
18-20 Miles Street
Oldham, Lancashire
OL1 3NU

1. Ideal for spore and fossil analysis.

2. Manufactured from clear acrylic in a 2-part design with a precision ground clearance to allow a variety of mesh grades to be tensioned, also incorporating a pouring lip for the convenience of decanting surplus liquid.

3. Made with a 75mm diameter sieving area and a wall thickness sufficient to provide a robust product that is not easily damaged.

4. Particularly suited to fine-end classification using synthetic meshes from 5-250 microns in aperture size.

5. Interchangeable meshes are available in the form of heat-sealed squares, please ask for our separate leaflet MSIVMESH.xls.



N O R M E S H

MiniSIV type A holder suitable for grades from 5-100 microns £61.00 each.

MiniSIV type C holder suitable for grades from 100-250 microns £61.00 each.

N.B. Prices subject to Post & Packaging and VAT (where applicable).

© 1998-2003 Mark Brown

